

QX laser is effective in clearing tattoos



Tattoo removal after first treatment and after fourth treatment

TATTOOS WERE ONE of the biggest fashion trends in the past decade. Now that many regret having one, demand for tattoo removal

is increasing. Tattoo removal by laser is the only effective and patient-friendly method that leaves no scarring.

Lasers deliver high peak power in short pulses, which shatter the ink particles or pigments without affecting the skin surface.

Fotona offers the Fotona QX, a Q-switched Nd:YAG and frequency doubled KTP Nd:YAG combo laser system, to remove even multi-coloured tattoos. Unlike comparable devices, the

Fotona QX is equipped with a zoom focus handpiece for both laser sources, eliminating the need to switch between handpieces during procedures.

The QX's high-energy capabilities and homogeneous laser beam profile allow larger laser spots, resulting in faster and safer treatments.

Website: www.fotona.si

Novatox tackles anti-ageing

A HIGH DOSAGE anti-ageing

