

Fotona

medical laser systems



	SP DYNAMIS Er:YAG+Nd:YAG	SP SPECTRO Er:YAG+Nd:YAG	XS DYNAMIS Er:YAG	XP DYNAMIS Nd:YAG
Laser Type	Er:YAG	Er:YAG	Er:YAG	
Wavelength	2940 nm	2940 nm	2940 nm	
Modalities	Full-field, Fractional V-SMOOTH, TURBO	Full-field, Fractional V-SMOOTH, TURBO	Full-field, Fractional V-SMOOTH, TURBO	
Energy / Power	3 J	3 J	3 J	
Scanner	F-Runner, S-Runner	F-Runner, S-Runner	F-Runner, S-Runner	
Laser Type	Nd:YAG	Nd:YAG		Nd:YAG
Wavelength	1064 nm	1064 nm		1064 nm
Modalities	LP,FRAC3,QCW	LP,FRAC3 QCW - optional		LP,FRAC3,QCW
Energy / Power	80 W	35 W		80 W
Scanner	S-11	S-11		S-11

Fotona
Chicago perfection

MULTI Well

ETG
Ergonomics
Technology
Group

JetPeel™
Fractional Laser

sanomed
medikal teknoloji

Bağdat Cad. Hasan Amir Sk. Dursoy İş Merkezi
No:4 Kat:3 Kızıltoprak / Kadıköy - İstanbul
Tel:0216 338 49 10 (pbx) Fax:0216 338 49 26
info@sanomed.com.tr - www.sanomed.com.tr